



InnovaSystems, Inc.

Technology Engineering

- **Automation to Support Metered Dose Inhaler Canister Cutting**
- Automates Canister Cutting for Total Can Content Analysis
- Accommodates Most Canister Sizes
- Reduces Injuries

MDI AC² – Automated Canister Cutter

The MDI AC²...

Automates the routine cutting of caps or containers in preparation for total container content and other tests.



Station accommodates a range of container sizes:



Is Easy to Use



6 Easy Steps:

1. Place depressurized canister in adapter
2. Lock canister in place
3. Close cover
4. Engage blade
5. Raise cover
6. Remove cut canister

MDI AC² includes

- MDI AC² Automated Canister Cutter Station
- One (1) Canister Adapter
- Users' Guide

Specifications

Space Requirements	Height	Width	Depth
Station	8.5in	7.4in	7.0in
	Power		
Utilities	110/220 VAC 50/60 Hz		

For more information contact

InnovaSystems, Inc.
1245 N. Church Street, Unit 6
Moorestown, NJ 08057
856-722-0410

www.InnovaSystemsInc.com
info@InnovaSystemsInc.com